


Client report

	No of participants	No of responded participants	Response percentage
Borrelia burgdorferi, antibodies, European origin, January, 1-2019	69	67	97.1 %

Summary

Overall success rate by samples




Summary	Own score	Max score	Own success rate	Difference	AVR success rate
Sample S001	8	8	100 %	1 %	99 %
Sample S002	8	8	100 %	0.7 %	99.3 %
Average:			100 %	0.9 %	99.1 %

History	Test nr.	Own success rate	Difference	AVR success rate
History not found				


Sample S001 |

Sample S001 success rate



Sample S001 results	Responded	Own score	Max score	Own success rate	Difference	AVR success rate	Count
	Clinical interpretation	4	4	100 %	0 %	100 %	78
	BorrAb	-	-	-	-	100 %	5
	BorrAbG	2	2	100 %	1.3 %	98.7 %	76
	BorrAbM	2	2	100 %	3 %	97 %	67
	Total:	8	8	100 %	1 %	99 %	226


Sample S001 Clinical interpretation



Sign of past or on-going infection
 Positive, referred for confirmation
 Our laboratory does not give a clinical interpretation to the physician


Clinical interpretation	Interpretation	Interpretation count	Own score	Max score	Own success rate	Difference	AVR success rate
	<input checked="" type="radio"/> Sign of past or on-going infection	30	4	4	100 %	0 %	100 %
	Positive, referred for confirmation	39					100 %
	Our laboratory does not give a clinical interpretation to the physician	9					-
	Total:	78	4	4	100 %	0 %	100 %

Sample S001 BorrAb




Positive

Sample S001 BorrAbG



Borderline Positive

Sample S001 BorrAbM



Negative Borderline Positive

BorrAb	Interpretation	Method	Interpretation count	Method count	Own score	Max score	Own success rate	Difference	AVR success rate
	Positive		5						100 %
		Immunitics C6-LYME-ELISA		5					
	Total:		5		-	-	-	-	100 %

OWN DEVICE: MET MANUALNA

BorrAbG	Interpretation	Method	Interpretation count	Method count	Own score	Max score	Own success rate	Difference	AVR success rate
	Borderline		1						0 %
		In-house IFA		1					
	<input checked="" type="radio"/> Positive		75		2	2	100 %	0 %	100 %
		Aesku Diagn AESKULISA Borrelia IgG ELISA		2					
		Aesku Diagn. Aeskublots Borrelia IgG		1					
		Biomedica Borrelia IgG ELISA Recomb.		1					
		bioMérieux VIDAS Lyme IgG		18					
		Diasorin LIAISON Borrelia IgG		24					
		Euroimmun Anti-Borrelia plus VlsE Elisa IgG		2					
		Euroimmun Anti-Borrelia Select ELISA IgG		1					
		Euroimmun Lyme screen ELISA IgG		1					
		Euroimmun Westernblot/ EUROLINE-WB IgG		3					
		In-house Borrelia IgG ELISA		1					
		Mikrogen ELISA Recomwell IgG/IgM		1					
		Mikrogen Line Blot Recomblot IgG		2					
		Mikrogen recomBead Borrelia IgG		1					
		<input checked="" type="radio"/> Mikrogen recomLine Borrelia IgG		1					
		Orgentec Borrelia IgG Recomb.		6					
		Siemens Enzygnost Borreliosis IgG		6					
		Testline BLOT-LINE Borrelia HGA IgG		1					
		Thermo Scientific IDEIA Borrelia Burgdorferi IgG		1					
		Virotech Borrelia afzelii + VlsE IgG Europe ELISA		1					
		Virotech Borrelia Europe plus TpN17 LINE IgG		1					
	Total:		76		2	2	100 %	1.3 %	98.7 %

OWN DEVICE: MET MANUALNA


BorrAbM	Interpretation	Method	Interpretation count	Method count	Own score	Max score	Own success rate	Difference	AVR success rate
	Negative		2						0 %
		Diasorin LIAISON Borrelia IgM		1					
		In-house IFA		1					
	Borderline		5						100 %
		Diasorin LIAISON Borrelia IgM		4					

Copyright © Labquality Oy
 Labquality does not permit any reproduction for commercial purposes of any portion of the material subject to this copyright. Labquality prohibits any use of its name, or reference to Labquality EQA program, or material in this report in any advertising, brochures or other commercial publications. Labquality EQA data do not necessarily indicate the superiority of instruments, reagents, testing equipments or materials used by participating laboratories. Use of Labquality EQA data to suggest superiority or inferiority of equipments or materials may be deceptive and misleading. Proficiency test results are handled confidentially. Labquality will not issue any statements to third parties of the performance of laboratories in external quality assessment schemes unless otherwise agreed.

	In-house Borrelia IgM ELISA		1					
<input checked="" type="radio"/> Positive		60		2	2	100 %	0 %	100 %
	Aesku Diagn AESKULISA Borrelia IgM ELISA		1					
	Aesku Diagn. Aeskublots Borrelia IgM		2					
	Biomedica Borrelia IgM ELISA Recomb.		1					
	bioMérieux VIDAS Lyme IgM		16					
	Diasorin LIAISON Borrelia IgM		14					
	Euroimmun Anti-Borrelia Select ELISA IgM		2					
	Euroimmun Lyme screen ELISA IgM		2					
	Euroimmun Westernblot/ EUROLINE-WB IgM		3					
	Mikrogen Line Blot Recomblot IgM		2					
	Mikrogen recomBead Borrelia IgM		1					
	<input checked="" type="radio"/> Mikrogen recomLine Borrelia IgM		1					
	Orgentec Borrelia IgM Recomb.		6					
	Siemens Enzygnost Borreliosis IgM		6					
	Testline BLOT-LINE Borrelia HGA IgM		1					
	Virotech Borrelia afzelii IgM ELISA		1					
	Virotech Borrelia Europe LINE IgM Line Immunoblot		1					
Total:		67		2	2	100 %	3 %	97 %


Sample S002 |

Sample S002 success rate



Sample S002 results	Responded	Own score	Max score	Own success rate	Difference	AVR success rate	Count
	Clinical interpretation	4	4	100 %	0 %	100 %	78
	BorrAb	-	-	-	-	100 %	5
	BorrAbG	2	2	100 %	1.3 %	98.7 %	76
	BorrAbM	2	2	100 %	1.5 %	98.5 %	67
	Total:	8	8	100 %	0.7 %	99.3 %	226


Sample S002 Clinical interpretation



Sign of past or on-going infection
 Positive, referred for confirmation
 Our laboratory does not give a clinical interpretation to the physician


Clinical interpretation	Interpretation	Interpretation count	Own score	Max score	Own success rate	Difference	AVR success rate
	<input checked="" type="radio"/> Sign of past or on-going infection	29	4	4	100 %	0 %	100 %
	Positive, referred for confirmation	40					100 %
	Our laboratory does not give a clinical interpretation to the physician	9					-
	Total:	78	4	4	100 %	0 %	100 %

Sample S002 BorrAb




Positive

Sample S002 BorrAbG



Borderline Positive

Sample S002 BorrAbM



Negative Borderline Positive

BorrAb	Interpretation	Method	Interpretation count	Method count	Own score	Max score	Own success rate	Difference	AVR success rate
	Positive		5						100 %
		Immunitics C6-LYME-ELISA		5					
	Total:		5		-	-	-	-	100 %

OWN DEVICE: MET MANUALNA

BorrAbG	Interpretation	Method	Interpretation count	Method count	Own score	Max score	Own success rate	Difference	AVR success rate
	Borderline		1						0 %
		In-house IFA		1					
	<input checked="" type="radio"/> Positive		75		2	2	100 %	0 %	100 %
		Aesku Diagn AESKULISA Borrelia IgG ELISA		2					
		Aesku Diagn. Aeskublots Borrelia IgG		1					
		Biomedica Borrelia IgG ELISA Recomb.		1					
		bioMérieux VIDAS Lyme IgG		18					
		Diasorin LIAISON Borrelia IgG		24					
		Euroimmun Anti-Borrelia plus VlsE Elisa IgG		2					
		Euroimmun Anti-Borrelia Select ELISA IgG		1					
		Euroimmun Lyme screen ELISA IgG		1					
		Euroimmun Westernblot/ EUROLINE-WB IgG		3					
		In-house Borrelia IgG ELISA		1					
		Mikrogen ELISA Recomwell IgG/IgM		1					
		Mikrogen Line Blot Recomblot IgG		2					
		Mikrogen recomBead Borrelia IgG		1					
		<input checked="" type="radio"/> Mikrogen recomLine Borrelia IgG		1					
		Orgentec Borrelia IgG Recomb.		6					
		Siemens Enzygnost Borreliosis IgG		6					
		Testline BLOT-LINE Borrelia HGA IgG		1					
		Thermo Scientific IDEIA Borrelia Burgdorferi IgG		1					
		Virotech Borrelia afzelii + VlsE IgG Europe ELISA		1					
		Virotech Borrelia Europe plus TpN17 LINE IgG		1					
	Total:		76		2	2	100 %	1.3 %	98.7 %

OWN DEVICE: MET MANUALNA

BorrAbM	Interpretation	Method	Interpretation count	Method count	Own score	Max score	Own success rate	Difference	AVR success rate
	Negative		1						0 %
		In-house IFA		1					
	Borderline		4						100 %
		Diasorin LIAISON Borrelia IgM		3					
		In-house Borrelia IgM ELISA		1					

Copyright © Labquality Oy
 Labquality does not permit any reproduction for commercial purposes of any portion of the material subject to this copyright. Labquality prohibits any use of its name, or reference to Labquality EQA program, or material in this report in any advertising, brochures or other commercial publications. Labquality EQA data do not necessarily indicate the superiority of instruments, reagents, testing equipments or materials used by participating laboratories. Use of Labquality EQA data to suggest superiority or inferiority of equipments or materials may be deceptive and misleading. Proficiency test results are handled confidentially. Labquality will not issue any statements to third parties of the performance of laboratories in external quality assessment schemes unless otherwise agreed.

<input checked="" type="radio"/> Positive		62		2	2	100 %	0 %	100 %	
	Aesku Diagn AESKULISA Borrelia IgM ELISA		1						
	Aesku Diagn. Aeskublots Borrelia IgM		2						
	Biomedica Borrelia IgM ELISA Recomb.		1						
	bioMérieux VIDAS Lyme IgM		16						
	Diasorin LIAISON Borrelia IgM		16						
	Euroimmun Anti-Borrelia Select ELISA IgM		2						
	Euroimmun Lyme screen ELISA IgM		2						
	Euroimmun Westernblot/ EUROLINE-WB IgM		3						
	Mikrogen Line Blot Recomblot IgM		2						
	Mikrogen recomBead Borrelia IgM		1						
	<input checked="" type="radio"/> Mikrogen recomLine Borrelia IgM		1						
	Orgentec Borrelia IgM Recomb.		6						
	Siemens Enzygnost Borreliosis IgM		6						
	Testline BLOT-LINE Borrelia HGA IgM		1						
	Virotech Borrelia afzelii IgM ELISA		1						
	Virotech Borrelia Europe LINE IgM Line Immunoblot		1						
	Total:		67		2	2	100 %	1.5 %	98.5 %